



™ Simple
Connected
Powerful



MTR 1 TB
531 hrs



UsbStorage
17hrs

10:16:35
2024/10/12




MTR

Touch Screen Video Recorder




Medical Imaging Made Easy™

Use Cases




Doctor's Office

Enabling seamless consultations or operation follow-ups with USB camera (PTZ) and audio kit for comprehensive documentation.



Microscopy

Makes it easy to record microscopic surgeries without obstructing the surgeon's view of the procedure.




USB Endoscopes ENT, Dental & Veterinary

Its an ideal solution for high-quality video capture in ENT, Dental & Veterinary, for comprehensive documentation.

MTR Specifications²

Device Type	Touch Screen Video Recorder	
Case	Aluminum Alloy	
Recording Resolution	3840x2160 ³ , 1920x1080, 1280x1024, 1280x720, 1024x768	
Video Recording Format	MP4: H.264 / HEVC H.265 High Efficiency Video Coding	
Image Format	JPG, PNG	
Display Dimensions	MTR133: 13.3" Diagonal	MTR156: 15.6" Diagonal
Touchscreen Type	Capacitive, 10-Point Multi-Touch	
Display Resolution	1920 x 1080 pixel	
Audio Input (recording)	MIC / USB Microphone / LINE / HDMI	
Audio Output	Built-In Loudspeakers / LINE / Headset / HDMI	
Internal Storage	1TB (SSD). Larger capacity can be ordered	
External Storage	USB Drives, Format FAT32 / NTFS / exFAT	
LAN	RJ-45, Ethernet, max 2.5 GB/s	
Dimensions (W x D x H)	MTR133: 306,6 x 184 x 32mm (12.8 x 7.2 x 1.3")	MTR156: 357,1 x 214,2 x 33mm (14.1 x 8.4 x 1.3")
Weight	MTR133: 1,1 kg (2.43 lbs)	MTR156: 1,3 kg (2.9 lbs)
Power	12V 3A external Power Supply	
Compliance and Certification	EN 60601-1 - Medical Equipment Safety EN 60601-1-2 - EMC, EN IEC 63000:2018 - Restriction of hazardous substances ICES003 (Canada), CE (MDR 2017/745/EU, WEEE, RoHS, FCC	
Order Information	MTR133 - Touch Screen Video Recorder (13.3") MTR156 - Touch Screen Video Recorder (15.6") AK4K - 4K Recording activation AKDICOM - DICOM software activation	

Optional add-ons Pay for what you really need



MediCapture MVR REMOTE APP



¹ Upgrade package option.

² These specifications are subject to change without notice.

³ For recording 4K input signals in 4K resolution, "4K recording" needs to be activated. Otherwise 4K input signals will be recorded in HD.

⁴ Will be available in the future.

MediCapture, Inc. - North America

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462, USA

+1 888 922 7887
www.medicapture.com
sales@medicapture.com

MediCapture Europe B.V.

Stationsplein 5 C2
6131 AT Sittard, The Netherlands

+31 46 208 1120
www.medicapture.com
sales.europe@medicapture.com

MediCapture APAC Inc.

3FL-6, No. 432, Section 1, Keelung Road,
Xinyi District 110, Taipei City, Taiwan, Republic Of China

+886 2 2345 3668
www.medicapture.com
sales.asiapacific@medicapture.com